

# HUGE

## PERFIT Tray Adhesive

Convenient and efficient impression accessory

## Tray Adhesive for all VPS Materials.


PERFIT Tray Adhesive guarantees excellent adhesion between VPS-based impression materials and metal or plastic trays. The product is remarkably easy to apply and straightforward to use. Impressions can be placed in the tray just 30 minutes after a brief application period of 3 minutes (max) .



# HUGE PERFIT Tray Adhesive

Indication :

To offer reliable adhesion between impression tray and silicone impression material.



● Suitable for A Silicone and C Silicone

● Easy to apply and handle



01




02

Short setting time after application

Drying Time < 3 minutes

● Does not create thickness

● Blue colour visual control



## Technical Parameters

PERFIT Tray Adhesive	
Drying Time	3:00
Bonding Strength (30min)	0.2N/mm <sup>2</sup>
Color	Blue
Shelf Life	3 Years

\*Data from HUGE Lab

## Order Information

PERFIT Tray Adhesive	
Standard	10ml Bottle



Shandong Huge Dental Material Corporation

Add /No.68 Shanhai Road, Donggang District, Rizhao City, Shandong Province, 276800, P.R. China.

Tel/ +86 (633) 227268 marketing@hugedental.com

www.hugedental.com



Facebook Huge Dental

Instagram Huge Dental

Youtube Huge Dental