

VinciSmile®

# HugeBond Universal FliPro

**8th generation**, one component self-etch light cure dental adhesive



## ► Universal Application

- Suitable for various substrates
- Capable of total-etch, self-etch and selective-etch
- Works with cement of light-cured, self-cured and dual cured
- Contains **MDP monomer**, delivering self-etching properties and higher bonding strength

## ► Gentle Experience

- Virtually no postoperative sensitivity
- Low occurrence of microleakage prevents secondary caries
- High moisture tolerance, providing consistently high bond strength on both dry and moist dentin

*Upgrade to flip-open cap  
Capacity increasing*

**5ML → 6ML**

## ► Reliable Bonding

- Bonding strength no lower than 20 MPa
- Low film thickness of 3µm
- Effectively seals over 90% dentin tubules
- Ultra-mild pH(PH=2.7) provides excellent dentin bond strengths and enamel bonds



# Versatile, Simple, Efficient...



- The latest innovation, **8<sup>th</sup> generation** adhesive.
- One-bottle technique, streamlines the workflow in practice.
- Featuring universality and versatility, compatible with multiple substrates in restorations.
- Plays a core role in clinical restorative solutions:
  - TopCEM Indirect Restorative Solution
  - TrusFIL Direct Restorative Solution
  - TopCEM Veneer Esthetics Solution




**27** CLINICAL EVALUATORS  
**660** TOTAL USES  
**86%** CLINICAL RATING




"EASY TO USE WITH MANY INDICATIONS."

## Dental Advisor Recommended:

Evaluation Summary:



Compared to Competitive Products:



VINCISMILE GROUP LLC

Add: 20524 Carrey Rd. Walnut, CA 91789, USA  
 Email: [marketing@vincismile.com](mailto:marketing@vincismile.com) Tel: +1 626 283 5808  
 Web: [www.vincismile.com](http://www.vincismile.com) Shop with us online at [store.vincismile.com](http://store.vincismile.com)

Follow @VinciSmile